

ROWHEELS® **REV** Order Form

| PROVIDER INFORMATION | | RTS INFORMATION | |
|-----------------------------|-------------|----------------------------|-------------|
| Provider: | | Name: | |
| Provider POC: | | Email: | |
| Account #: | | Phone #: | |
| Email: | | Address: | |
| Phone #: | | City: | State: ZIP: |
| Fax #: | | Ship to information | |
| Address: | | Company (if applicable): | |
| City: | State: ZIP: | Name: | |
| Client: | | Phone #: | |
| Rowheels information | | Address: | |
| Ship date: | | City: | State: ZIP: |
| Order #: | | Ship via: | |
| User Information | | Wheelchair Information | |
| Weight: | Height: | Wheelchair Make: | Model: |
| Med. Diagnosis: | | Camber Angle: | Other: |

Rowheels are covered through HCPSC billing code E2227

Please email or fax to Rowheels customer service at orders@rowheels.com or (608) 719-0010



Key 1: MODEL

| MODEL | KEY CODE | KEY INPUT | |
|------------------------|----------|-----------|-----------------|
| REV-LX (Low gear)..... | RWLX | | \$5,400 (E2227) |
| REV-HX(High gear)..... | RWHX | | \$5,400 (E2227) |



The Right Way to Roll.™

www.rowheels.com

2895 Commerce Park Drive
Fitchburg, WI 53719

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

ROWHEELS® **REV** Order Form

Key 2: COLOR SCHEME

| COLOR SCHEME | KEY CODE | KEY INPUT |
|----------------------------------|--|-----------|
| Red/Black..... | 0.....STD | |
| Black/Black..... | 1.....STD | |
| White/Black..... | 2.....STD | |
| Red, Black or White/Custom | C.....RFQ (Custom hydrographic on tri-adapter) | |



Red/Black



Black/Black



White/Black

Key 3: TORQUE ARM

| TORQUE ARM | KEY CODE | KEY INPUT |
|---------------------|------------|-----------|
| 8" Torque Arm..... | 8.....STD | |
| 10" Torque Arm..... | 10.....STD | |

Key 4: HANDRIM

| HANDRIM | KEY CODE | KEY INPUT |
|-------------------------------|---|-----------|
| Out-Front Natural-Fit..... | N.....\$335 | |
| Out-Front Natural-Fit LT..... | NL.....\$365 | |
| Out-Front SURGE..... | S.....\$425 | |
| Out-Front SURGE LT..... | SL.....\$475 | |
| Q-Grip..... | Q.....\$315 (not available on 22" or 26" wheels) | |
| RibGrips..... | G.....\$260 (only available on 24" wheels & in blue, pink, black and red) | |
| Projection-8 (45°)..... | P.....\$375 | |



The **Natural-Fit** has an oval profile that increases gripping surface area by 80%. The optional thumb grip component is not available for use on Rowheels.

The **Surge** has the same oval profile along with a rubber strip along the top surface of the oval for increased friction while leaving the side surface smooth for braking.

Q-Grips are ideal for anyone with limited hand function or anyone that prefers the feel & functionality of a fully-coated handrim. Research shows the Q-Grip to be 59.1% more grippy than a vinyl handrim.

RibGrips are ideal for those looking to push their wheelchair with more ease and comfort. The soft rib discs benefit users with poor hand/finger grip.



The Right Way to Roll.™

www.rowheels.com

2895 Commerce Park Drive
Fitchburg, WI 53719

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

ROWHEELS® **REV** Order Form

Key 5: AXLE LENGTH

| AXLE LENGTH | KEY CODE | KEY INPUT |
|----------------------|-----------|--|
| Length 0 (2.0")..... | 0.....STD | <i>*Length 0 fits most Tilite, Kimobility, Pride frames</i> |
| Length 1 (2.1")..... | 1.....STD | <i>*Length 1 fits most Colours, Motion Composites frames</i> |
| Length 2 (2.2")..... | 2.....STD | <i>*Length 2 fits most Quickie frames</i> |
| Length 3 (2.3")..... | 3.....STD | <i>*Length 3 fits Top End Crossfire/T6 frames</i> |
| Length 4 (2.4")..... | 4.....STD | <i>*Length 4 fits Invacare ProSpin/Patriot Plus, Top End All-Terrain</i> |
| Length 5 (2.5")..... | 5.....STD | |

Key 6: TIRE

| TIRES | KEY CODE | KEY INPUT |
|---|--|-----------|
| No tire..... | N.....NA | |
| Schwalbe RightRun (Black/Grey)...R..... | STD | |
| Schwalbe Marathon Plus (Black)...M..... | \$130 | |
| Shox Solid (Black).....S..... | \$59 | |
| Kenda Nevegal.....X..... | \$78 (not available on 22" & 26" wheels) | |



The Schwalbe **RightRun** tires incorporate their most effective SmartGuard puncture protection layer, hand-friendly 2Grip side walls, and attractive two color compound. Robust Black'n Roll tread center compound that is non-marking



The Schwalbe **Marathon Plus** is the benchmark for puncture-proof pneumatic tires on the wheelchair. The all-new version comes with improved rolling resistance along with hand-friendly 2Grip side walls, and the robust Black'n Roll tread center compound.



SHOX Solid Wheelchair Tires are airless so you don't have to worry about punctures and are a firm 110psi equivalent.



The Kenda **Nevegal** is the optimal tire for an assortment of ever-changing terrain. It is highly adaptable and excels under various conditions, from loose rock, tall grass, muddy roots, single track and high speed descents, and many more.

Key 7: WHEEL SIZE

| WHEEL SIZE | KEY CODE | KEY INPUT |
|------------------|-----------|-----------|
| 22" (489mm)..... | 2.....STD | |
| 24" (540mm)..... | 4.....STD | |
| 25" (559mm)..... | 5.....STD | |
| 26" (590mm)..... | 6.....STD | |



The Right Way to Roll.™

www.rowheels.com

2895 Commerce Park Drive
Fitchburg, WI 53719


PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com

ROWHEELS® **REV** Order Form

KEY 8: SPOKE GUARD

| SPOKE GUARD | KEY CODE | KEY INPUT |
|---------------|----------|-----------|
| Carbon Fiber* | S (STD) | |
| Clear | C (STD) | |
| Custom | Z (RFQ) | |

*Guard is covered in a vinyl carbon fiber wrap, Not actual carbon fiber



FINAL ORDER KEY CODE

Enter, in order, the seven option keys selected on the previous page here:

KEY#1 KEY#2 KEY#3 KEY#4 KEY#5 KEY#6 KEY#7 KEY#8

Sample Final Order Key code - **RWLX-0-8-N-1-R-4**


Description - REV 1, Red/Black, 8" Torque Arm, Natural-Fit, Length 1, Right Run Tire, 24"

REV¹ WHEELCHAIR ADAPTER INFO

In order for REV wheels to function properly, our universal frame clamp latches must be installed on the wheelchair's frame. If the chair to be used has a tubular frame with a diameter between **0.75" (19.05mm) and 1.375" (34.93mm) or is a Quickie Q7**, select "Compatible Frame", otherwise select "Not Compatible" and provide us with the wheelchair make, model and frame dimensions in the "Notes" section below. A decision regarding whether a custom adapter can be designed will be made. An RFQ can then be submitted for the cost of the custom adapter.

Compatible Frame

Not Compatible



NOTES:



The Right Way to Roll.™

www.rowheels.com

2895 Commerce Park Drive

Fitchburg, WI 53719

PHONE: (608) 268-9670 | FAX: (608) 719-0010 | support@rowheels.com